

# Pregnancy Test Early, Docalon 2 Count Clear and Accurate Results Over 99% Accurate HCG Early Pregnancy Test Detection Individually Sealed Review-2021



[Read Customer Reviews](#)

Detects all forms of the pregnancy hormone commonly found in urine. Suggestion: Prepare a cup, the height of the urine in the cup is about 2.3 cm approximately. (You can measure it from the outside of the cup, the approximate height is when the urine reaches the fifth mark, which is the highest mark). There should not be too much of urine in the cup. If there is too much, it will dilute the "gold" in the product and the display will have no lines. It can't be too little, because if there is too little urine, the product cannot absorb the urine and react with "gold". The absorption time of the product in the cup is about 30-50 seconds (if you don't count the seconds, you can wait for the display to appear pink, and then you can take it out after 2 more seconds). Keep it still and wait for the result, usually it will appear within 1 minute. Suggestion: If you don't understand, you can also watch our video, which has a detailed explanation of the process. Use any time of day  
Two pink lines pregnant, one pink line not pregnant. A pink and a pale red are also pregnant  
Each box contains: two bags and desiccant in each box. Store at 36°-86°F (2°-30°C). Do not freeze. For in vitro diagnostic use. Not for internal use. Keep out of the reach of children. Read package insert before use

That our pregnancy test sticks offer you the best value and quality in the market. If you have any questions, please feel free to contact the seller. We will reply to you within 12 hours. Our target is to offer accurate test to make you get accurate result. Early knowing will make you pay attention to pregnancy in advance and better protect your pregnancy. If your test is positive, please consult your doctor.

Warm reminder Common drugs (such as cold medicine, antibiotics, painkillers, etc.) will not affect the accuracy of the test. However, injecting or taking HCG hormones and other AIDS drugs can affect pregnancy. Certain conditions and diseases that can disrupt hormone levels can affect the menstrual cycle, leading to different hormone readings. Your results can be affected by menopause, taking hormones or steroids, birth control pills, polycystic ovary syndrome, hyperthyroidism, and endocrine disorders.

Docalon HCG Test Kits (2 individual package).

Usually only 1 test. If you doubt, please use morning urine to retest and confirm. Please trust once you get the result.

It is accurate just like urine test in Hospital.

The precision rate exceeds 99% when test the typical HCG level by this testing kit. Why is it 99%? As the body status and condition vary from people to people. Some cases will not reach the effect. but that is only individual person. When the urine HCG reaches 25mlu/ml, it shows positive. The HCG for person without pregnancy is 0--5mlu/ml, it is negative.

Detects four Major Forms of hCG.

Docalon early HCG test kit determine all common pregnant hormone urine. The urine sample must reach the upper hole so that the sample can be fully absorbed to get more accurate result.

Reliable test is a must.

Nothing is more exciting than that of welcoming the coming of new baby. Our Docalon with 99% precision will let you share your joy and happiness with your relatives early. Review 2021, feedback 2021, promo code, discount code 2021, buy, picture, description, sale, price comparison, cheap, cheapest, value for money.